

MagPure Blood RNA Kit C

Introduction

This product supplies a simple and rapid extraction of total RNA including microRNA from Blood, buffy coat, bone marrow, Cell suspension and other body fluids. The kit is based on superparamagnetic particles purification technology, no phenol-chloroform extraction or alcohol precipitation. Purified RNA is ready for downstream applications such as RT-PCR, virus RNA testing and so on.

Principle

The Kit can be used for both manual extraction process and automatic nucleic acid extraction machines. This Kits is suitable for extracting RNA from $\leq 1 \times 10^6$ cells suspension, 50 μ l Whole Blood, 50 μ l buffy coat, 20 μ l bone marrow. This product is based on the purification method of high binding magnetic particles. The sample is lysed and digested under the action of lysate and Protease. RNA/DNA is released into the lysate. After adding magnetic particles and binding solution, DNA/RNA will be adsorbed on the surface of magnetic particles, and impurities such as proteins will be removed without adsorption. The adsorbed particles were digested using DNase and washed with washing buffer to remove proteins and impurities, washed with ethanol to remove salts, and finally RNA was eluted by Elution Buffer.

Kit Contents

Cat.No.	R661101C	R661102C	R661103C
Purification times	48 Preps	96 Preps	480 Preps
MagPure RNA Particles	1.0 ml	1.5 ml	6 ml
Proteinase K	24 mg	48 mg	220 mg
Protease Dissolve Buffer	1.8 ml	1.8 ml	15 ml
DNase I	600 μ l	2 x 600 μ l	10 x 600 μ l
DNase Buffer	15 ml	15 ml	60 ml
Buffer AL	10 ml	10 ml	60 ml
Buffer MCB*	9 ml	9 ml	30 ml
Buffer MW1*	13 ml	22 ml	110 ml
Buffer MW2*	6 ml	20 ml	100 ml
RNase Free Water	10 ml	15 ml	120 ml

Storage and Stability

MagPure RNA Particles and Proteinase K should be stored at 2–8°C upon arrival. DNase I should be stored at -20°C. However, short-term storage (DNase I up to 1 weeks, MagPure RNA Particles and Proteinase K up to 8 weeks) at room temperature (15–25°C) does not affect their performance. The remaining kit components can be stored at room temperature (15–25°C) and are stable for at least 18 months under these conditions.

Materials and Equipment to be Supplied by User

- 100% ethanol
- Dilute Buffer MW1 with 17ml (48 Preps), 28ml (96 Preps) or 140ml (480 Preps) 100% ethanol and store at room temperature
- Dilute Buffer MW2 with 24ml (48 Preps), 80ml (96 Preps) or 400ml (480 Preps) 100% ethanol and store at room temperature
- Dilute Buffer MCB with 21ml (48 Preps), 21ml (96 Preps) or 70ml (480 Preps) isopropanol and store at room temperature
- Dissolve the Proteinase K (40mg/ml) with 0.6ml (48 Preps), 1.2ml (96 Preps) or 1.1ml (480 Preps) protease Dissolve Buffer to the Proteinase K and store at -20~8°C.

Manual or Liquid station protocol

1. Pipet 10µl Proteinase K(40mg/ml) and 10µl MagPure RNA Particles into the bottom of a 1.2ml well Plate (2.2ml).
2. Add 50µl sample to the 96 well plate.
3. Add 70µl Buffer AL to the sample and mix by shaking at 900~1000rpm for 15 min.
4. Add 140µl Isopropanol to the sample and mix by shaking for 5min. **Place the deep well plate on an Magnet Plate and allow beads to separate for 2 minutes.** With the plate on the Magnet Plate, perform the aspiration, and then discard the supernatant from the plate.
5. **Add 200µl Buffer MW1 and shaking for 2 minute to re-suspend the particles.** Place the tube to the magnetic rack for 1 minute, then remove the supernatant.
6. Leave the plate on the magnetic separation device. Wait 1 minute. Remove residual liquid with a pipettor. Dry the MagPure RNA Particles for an additional 3 minutes.

7. **Add 70µl DNase Mixture (40µl DNase Buffer + 10µl DNase I) to the sample.** Mix by shaking at 600-900rpm for 10~15 minutes.
DNase I and DNase Buffer can be premixed.
8. **Add 140µl Buffer MCB to the sample, shaking at 900-1200rpm for 6 minutes.** Place the tube to the magnetic rack for 1 minutes, then remove the supernatant.
9. **Add 200µl Buffer MW1 and shaking 900~1200rpm for 1 minute to re-suspend the particles.** Place the tube on the magnetic rack for 1 minutes, then remove the supernatant.
10. **Add 200µl Buffer MW2 and shaking 900~1200rpm for 1 minute to re-suspend the particles.** Place the tube on the magnetic rack for 1 minutes, then remove the supernatant.
11. Repeat step 10 once.
12. Leave the plate on the magnetic separation device. Wait 1 minute and remove residual liquid with a pipettor.
13. Dry the Mag-Pure Particles for an additional 10 minutes.
14. **Add 50µl RNase Free Water to sample and mix by shaking for 5 minutes.** Place the tube to the magnetic rack for 3 minutes.
- 15.** Transfer the supernatant containing the purified RNA to a new Plate and store RNA at -80°C.

KingFisher or similar Extractor isolation:

1. Add the Reagents/sample to the well of f the deep well plate according to the table below.

Name of the Plate	Pre-loaded reagents	Addition before use
Sample plate	1. Pipet 10µl Proteinase K and 10µl MagPure RNA Particles 2. Pipet 50µl of the sample into the well of plate. 3. Pipet 70µl Buffer AL. 4. Mix well and incubate for 15 mintues 5. Add 140µl Isoprpoanol	
Wash Plate 1	200µl Buffer MW1, Put in 96 magnetic Tip	
DNase	60µl DNase Buffer and 10µl DNase I	
Wash Plate 2	200µl Buffer MW2	
Wash Plate 3	200µl Buffer MW2	
Elution plate	50µl RNase Free Water	

2. Place a 96 tip comb for deep well magnets on Wash Plate 1.
3. Start the R6611C_Flex protocol with the KingFisher Flex 96 and load the plates.
4. **Add 140µl Buffer MCB to the DNase plate during the dispense step.**
5. Place the DNase plate back into the instrument and press Start. After the pause, the protocol will continue to the end.
6. After the run is completed, remove the plates and store the purified total RNA.

Troubleshooting Guide

1. Low RNA yields

- **Incomplete resuspension of MagPure RNA Particles:** Resuspend the MagPure RNA Particles by vortexing before use.
- **Loss of MagPure RNA Particles during procedure:** Be careful not to remove the MagPure RNA Particles during the procedure.
- **MagPure RNA Particles not resuspended during binding:** Vortex vigorously for 2 minutes after addition of Buffer MCB.
- **Eluate contains residual ethanol:** Ensure that the wash flow-through is drained from the collection tube and that the column is then centrifuged at $>12,000 \times g$ for 1 min.

2. Low A260/A280 value

- **Water used to dilute RNA for A260/A280 measurement:** Use 10 mM Tris-Cl, pH 7.5, not RNase-free water, to dilute the sample before measuring purity..