

Product Name】 MagPure Universal RNA Precast Kit (Auto Pure 32)

【Product Specification】 96 Preps/Kit

【Intended Use】

This product is suitable for rapid extraction of RNA from tissue, cells, blood, s and other clinical samples. RNA can be used directly for RT-PCR, quantitative RT-PCR and so on.

【Principle】

The Kit combines the speed and efficiency of silica-based technology with the convenient handling of magnetic particles for purification of total RNA. Samples are lysed and RNA is purified from lysates in one step through its binding to the silica surface of the particles in the presence of a chaotropic salt. The particles are separated from the lysates using a magnet and DNA is removed by treatment with RNase-free DNase I. The magnetic particles are efficiently washed, and RNA is eluted in RNase-free water.

【Main Composition】

Cat.No	Precast Reagent	IVD3020-TL-06	IVD3020-TL-06-00
Purification times		96 Preps	16 Preps
DNase I		2 x 600 μ l	600 μ l
DNase Buffer		60 ml	15 ml
RTL Lysis Buffer		60 ml	15 ml
Buffer MCB		18 ml	3 ml
96-Tip		12 PCS	2 PCS
2.0ml V-bottom plate	Row 1/7: 500 μ l Buffer MCB	6 Plates	1 Plate
	Row 2/8: 500 μ l Buffer MW1		
	Row 3/9: empty		
	Row 4/10: 30 μ l Magpure RNA Particles, 500 μ l Buffer MW2		
	Row 5/11: 900 μ l Buffer MW2		
	Row 6/12: 80 μ l RNase Free Water		

【Storage conditions and validity】

DNase I should be shipped with ice pack and stored at -20°C after arrival. However, short-term storage (up to 1 week) at room temperature (15–25°C) does not affect its performance. The remaining kit components can be stored at room temperature (15–25°C) for 18 months.

【Applicable Instrument】

Nucleic Acid Extraction Machine such as Auto Pure 32 (Allsheng) or similar.

【Preparation before Use】

- (Optional) 2-mercaptoethanol can be added to an aliquot of RTL Lysis Buffer before use. Add 20 μ l 2-mercaptoethanol per 1ml RTL Lysis Buffer. This mixture can be stored for 2 weeks at room temperature
- Add 10ml (16 Preps) or 42 ml (96 Preps) isopropanol to the bottle of MCB.

【Part 1: Sample Preparation】

- Cells Grown in Suspension (no more than 5×10^6 cells)**
Pellet cells by centrifugation. loosen the cell pellet thoroughly by flicking the tube. Add 500~600 μ l of RTL Lysis Buffer and vortex vigorously.
- Cells Grown in Monolayer (no more than 5×10^6 cells)**
Lyse cells directly in a culture dish by adding 600 μ l of RTL Lysis Buffer and passing the cell lysate several times through a blue pipette tip. Always use more RTL Lysis Buffer if in the lysate is too viscous to aspirate with a pipette.
- Animal Tissue (Do not use more than 20 mg Tissue)**
Homogenize no more than 20mg tissue samples in 600 μ l RTL Lysis Buffer using an appropriate mechanical homogenizer. Alternatively one can pulverize tissue in liquid nitrogen with mortar and pestle and transfer the powder to a clean 1.5 ml microcentrifuge tube and add RTL Lysis Buffer. After lysate, centrifuge at 14,000 x g for 3 minute at room temperature.
- Plant Tissue (Do not use more than 100 mg Tissue)**
Disruption Plant sample by liquid nitrogen and transfer up to 100mg power to 1.5ml Tube. Add 600 μ l RTL Lysis Buffer to the sample and mix well by vortexing vigorously. Centrifuge at 14,000 x g for 3 minute at room temperature.
- Whole Blood (Do not use more than 1.5ml blood)**
Separate the leukocytes cell from 0.5~1.5ml Whole blood. Resuspend the leukocytes pellets completely by 50 μ l Buffe PBS. Add 450 μ l of RTL Lysis Buffer and vortex vigorously to lyse.

- f. **Trizol/MagZol Reagent (without chloroform extraction):**
Following Trizol/MagZol Reagent protocol, add 600µl MagZol Reagent to lyse the sample. After lysate, centrifuge at 12,000 x g for 10 minutes at 2-8°C.
- g. **Trizol/MagZol Reagent (with chloroform extraction):**
Following Trizol/MagZol Reagent protocol, add 1ml Trizol/MagZol Reagent to lyse the sample. After lysate, add 200µl chloroform to the lysate and centrifuge at 12,000 x g for 15 minutes at 2-8°C.














【Part 2: Auto Pure 32 nucleic acid extractor operation】

1. Take out the required components of the kit.
2. Inverting the Plate several times to re-suspend the magnetic beads. Remove the sealing bag and sealing film.
3. **Add 400~450µl the lysate or the supernatant into each well of Raw 1/7.**
4. **Add 290µl DNase Buffer and 10µl DNase I into each well of Raw 3/9.**
DNase Buffer and DNase I can mix before adding in to the raw.
5. Load the plates and AS-tip on the machine.
6. Turn on the machine and start the IVD3020-TL-06 protocol.
7. The program pause at about 20 minutes. Take out the plate and **add 450µl Buffer MCB to each well of Raw 3/9.**
8. Place the plate back into the instrument and continue the program.
9. After the run is completed, remove the plates tips
10. Transfer the purified total RNA into new 1.5ml centrifuge tubes and store at -20 °C.

【Basic Information】

	Guangzhou Magen Biotechnology Co., Ltd.		
	Room 401, Building D, No. 7, Jingye 3rd Street, Yushu Industrial Park, Guangzhou Hi-Tech Industrial Development Zone, Huangpu District, Guangzhou, 510663, China www.magen-tec.com +86-20-3855 5004 info@magen-tec.com		

【Explanation of Marks】

	The product is used in vitro, please don't swallow		Please don't reuse it
	Validity		Please read the instruction book carefully before using
	Warning, please refer to the instructions in the annex		Manufacturer
	Temperature scope within which the product is reserve		Batch number
	European union authorization representativ		Keep dry
	Avoid overexposure to the sun		Don't use the product when the package is damaged
	The product meets the basic requirements of European in vitro diagnostic medical devices directive 98/79/EC		