

# Initial Library Preparation Bead Evaluation

## I. Objective

To compare the performance of three batches of MagPure A4 XP (Magen #XP-50/XP-500) with AMPure XP (Beckman #A63881).

## II. Materials

Comparison between three batches of MagPure A4 XP (Magen) and AMPure XP (A63881):

- AMPure XP (A63881)
- Three batches provided by Magen:
- HG050200, HE150200, HG070200

## III. Experimental Design

1. For each batch of library preparation beads, three libraries were constructed using:

- Single size selection method
- Double size selection method

2. A total of 32 cfDNA samples were pooled into one tube, mixed thoroughly by vortexing, and then aliquoted into 28 tubes, each containing 30  $\mu\text{L}$ .

3. Library yields were compared across different bead batches using the same library preparation method.

## IV. Results

1. cfDNA Concentration

0.208 ng/ $\mu\text{L}$

2. Library Construction Results

Bead Batch	Source	Library Prep Method	Sample ID	Adapter	Library Concentration (ng/μL)	Mean (ng/μL)
AmPure XP	Library preparation beads currently in use in the laboratory	Single size selection method	1	4	8.73	6.35
			2	5	6.24	
			3	6	4.08	
HG050200	Magen QC batch beads		4	7	8.23	7.99
			5	8	7.57	
			6	9	8.16	
HG070200			10	4	9.1	8.32
			11	5	6.77	
			12	4	9.1	
HE150200			13	4	9.1	8.32
			14	5	6.77	
			16	7	8.7	
AmPure XP	Library preparation beads currently in use in the laboratory		17	8	4.32	3.39
			18	9	3.05	
			19	10	2.81	
HG050200	Magen QC batch beads		20	11	3.31	3.15
		21	12	2.84		
		22	13	3.31		
HG070200		23	14	2.99	3.71	
		24	15	3.46		
		25	4	4.69		
HE150200		26	5	4.01	3.72	
		27	6	3.54		
		28	7	3.62		

# Microscale Recovery Efficiency Evaluation

## Marker Recovery Experiment

### Objective

To evaluate the recovery efficiency of XP magnetic beads at different ratios using a **25 ng, 50 bp purified DNA marker**.

### Experimental Procedure

1. Add **20 ng of purified marker** into a **1.5 mL microcentrifuge tube**.
2. Adjust the volume to **100  $\mu$ L** with nuclease-free water.
3. Add XP magnetic beads (different batches) at ratios of **0.7 $\times$ , 0.8 $\times$ , 0.9 $\times$ , 1.2 $\times$  and 1.5 $\times$** , respectively.
4. Vortex to mix thoroughly for **1 min**.
5. Incubate at room temperature for **10 min** to allow binding.
6. Place the tube on a magnetic stand for **5 min**.
7. Carefully remove and discard the supernatant.
8. Wash the beads twice with **500  $\mu$ L of 70% ethanol**.
9. Briefly centrifuge and remove any residual ethanol.
10. Air-dry the beads for **10 min**.
11. Elute DNA in **50  $\mu$ L of EB buffer**.

		qubit	Yield	Recovery
A	Original sample dilution	0.502	0.0251	100%
AMPure XP	0.7	0.0975	0.004875	19%
	0.8	0.205	0.01025	41%
	0.9	0.274	0.0137	55%
	1.2	0.477	0.02385	95%
	1.5	0.475	0.02375	95%
HG050200	0.7	0.112	0.0056	22%
	0.8	0.219	0.01095	44%
	0.9	0.284	0.0142	57%
	1.2	0.468	0.0234	93%
	1.5	0.486	0.0243	97%
HG070200	0.7	0.106	0.0053	21%
	0.8	0.216	0.0108	43%
	0.9	0.286	0.0143	57%
	1.2	0.475	0.02375	95%
	1.5	0.485	0.02425	97%