Nuclease Free Water (DNase/RNase Free)

Description

Nuclease Free Water is DEPC-free ultrapure water, suitable for experiments that are sensitive to DEPC residues.

Ordering Informatoin

CAT. No.	Product Name	Package
C454	Nuclease Free Water	500 ml

Product Parameters

Application	Dissolution of RNA
Packaging	Polypropylene plastic bottle
Composition	Ultrapure water
Package	500 ml
DNase	Not detected
RNase	Not detected
Sterilization	Autoclaved (121°C, 20 minutes)
Storage Conditions	Store at room temperature. For long term storage, store at 2–8°C
Shelf Life	12 months
Preparation	Prepared with ultrapure water